



INTENDED USE

Wool products can be used in the rehabilitation and prophylaxis of people affected by: rheumatoid and degenerative diseases of the upper spine, ankylosing muscle inflammation, the effects of injuries (fractures, dislocations, deformities, contractures, joint instability) and by people prone to excessive sweating, allergic to feathers and mites.

USE

Through the stabilization in the cervical region and the heating of this area, wool products eliminate the forced tension of the shoulder and cervical spine muscles during a long journey and allow the correct positioning of the head during lying and sleeping.

The products can also be used as supports of neck and nape during rehabilitation exercises and after the end of rehabilitation, in order to prolong the thermal effect and for everyday use.

CONTRAINDICATIONS

- Should not be used by people who are allergic to sheep's wool.
- Should not be used by people who have damaged skin or wounds.

WARNINGS

In case of babies, small children and people sensitive to wool, it is recommended to use woolen products coated in cotton fabric covers. It will protect the users against loose wool fibers getting to the eyes, nose and throat and the associated unpleasant irritating effect.

WASHING INSTRUCTIONS

- Hand wash at temperatures up to 30° C with the use of wool detergents containing lanolin.
- Do not bleach with chlorine-evolving compounds.
- Hang out to dry.
- Do not spin.
- Do not iron.
- Can be dry cleaned.



Rehabilitation products made of wool fibers may have a specific smell.
It is recommended to air the products until its intensity is reduced.